

Drug Of Abuse - URINA

Testul **One-Step Drug of Abuse** este un imunotest cromatografic cu flux lateral pentru detectarea calitativa a mai multor medicamente si metaboliti ai medicamentelor in urina umana la o concentratie limita.

AMP Amphetamine 500ng/ml; 1000ng/ml						
BAR Barbiturates 300ng/ml	BUP Buprenorphine 10ng/ml	BZO Benzodiazepines 200ng/ml ; 300ng/ml				
CAF Caffeine 10000ng/ml	COC Cocaine 100ng/ml ; 200ng/ml ; 300ng/ml	COD Codeine 300ng/ml	COT Cotinine 100ng/ml ; 200ng/ml ; 500ng/ml			
EDDP Methadone CLO Clonazepam Metabolite						
150ng/ml ; 400ng/ml	100ng/ml	1000ng/ml	200ng/ml			
KET Ketamine 500ng/ml ; 1000ng/ml	K2 Synthetic marijuana 50ng/ml ; 200ng/ml	LSD Lysergic Acid Diethylamise 20ng/ml ; 50ng/ml	6-MAM O6-Monoacetylmorphine 10ng/ml	MCAT Methcathinone 500ng/ml		
MET (mAMP) Methamphetamine 300ng/ml ; 500ng/ml ; 1000ng/ml	MDMA Ecstasy 500ng/ml	MDPV Methylenedioxy Pyrovalerone 1000ng/ml	MOP/TRA Tramadol / 100ng/ml ; 200ng/ml ; 300ng/ml	MPD Methylphenidate 300ng/ml ; 1000ng/ml	MTD/MQL Methadone/Methaqualone 300ng/ml	
OPI Opiate 300ng/ml ; 2000ng/ml	OXY Oxycodone 100ng/ml	PCP Phencyclidine 25ng/ml	PPX Propoxyphene 300ng/ml	TCA Tricyclic Antidepressants 1000ng/ml	THC Marijuana 25ng/ml; 50ng/ml	ZOL Zolpidem 50ng/ml



PROCEDURA TESTULUI SI INTERPRETAREA

Banda

Caseta

Dipcard

Recipient

3~5mins

Drug Of Abuse - Saliva

Testul **One-Step Drug of Abuse** este un imunotest cromatografic cu flux lateral pentru detectarea calitativa a mai multor medicamente si metaboliti ai medicamentelor in saliva umana la o concentratie limita.


Informatii Generale

- ◆ Numai pentru uz diagnostic in vitro
- ◆ Perioada de valabilitate: 24 luni
- ◆ A se depozita sigilat la 4~30°C (40~86°F)
- ◆ Proba: Saliva
- ◆ Cititi rezultatele 3~5 minutes

AMP Amphetamine 50ng/ml	BAR Barbiturates 300ng/ml	' BUP Buprenorphine 10ng/ml	BZO Benzodiazepines 50ng/ml
	coc Cocaine 20ng/ml	COT Cotinine 50ng/ml	KET Ketamine 100ng/ml
	MET (mAMP) Methamphetamine 50ng/ml	MOP Morphine 15ng/ml; 40ng/ml	MTD Methadone 75ng/ml
		OPI Opiate 40ng/ml	OXY Oxycodone 50ng/ml
		PCP Phencyclidine 10ng/ml	THC Marijuana 12ng/ml; 25ng/ml; 50ng/ml



PROCEDURA TESTULUI SI INTERPRETAREA



HANGZHOU CLONGENE BIOTECH CO., LTD.

Add: No.1 Yichuang Road, Yuhang Sub-district, Yuhang District, Hangzhou 311121, China

Tel: +86(0)57188262120 Fax: +86(0)57188261752

Web: www.clongene.com E-mail: marketing@clongene.com





杭州隆基生物技术有限公司
Hangzhou Clongene Biotech Co., Ltd.

Tel: +86-571-88262120
Fax: +86-571-88261752

Web: www.clongene.com
Email: marketing@clongene.com

DECLARATIE DE CONFORMITATE CE

Denumirea si adresa fabricii Hangzhou Clongene Biotech Co., Ltd.
No.1 Yichuang Road. Yuhang Sub-district
Yuhang District
311121 Hangzhou
China

Declaram pe propria noastra raspundere:

Dispozitivul medical: Multi-Drug 10 Test Dip Card Droguri-DDOA1a04

Din clasa: Altele
Conform articolului 9 al Directivei 98/79/EC

Indeplineste prevederile Directivei 98/79/CE si transpunerilor acesteia in legislatia nationala aplicabila.
Declaratia este valabila in legatura cu "raportul de inspectie final" al dispozitivului.

Procedura de evaluare a conformitatii: Directiva 98/79/EC Annex III, cu exceptia Sectiei 6

Standarde aplicabile: EN ISO 13485:2016 EN ISO 15223-1:2016
EN ISO 23640:2015 EN13612:2002/AC:2002 EN ISO
EN 13975:2003 14971:2012 EN ISO 18113-2:2011
EN ISO 18113-1:2011
EN 62366-1:2015

Numele si adresa reprezentantului autorizat: Shanghai International Holding Corporation GmbH (Europe) Eiffestrasse 80
20537 Hamburg Germania

Hangzhou. 9 August 2021

Locul, Data:



Numele si functia

Shi Jinqiang, Reprezentant Legal